

COVID-19 Testing Solutions

Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2)

SARS-CoV-2 IgM/IgG Antibody

Assay Kit by Colloidal Gold Method

20 Tests/Kit

Serum/Plasma: 10 μ L

Whole Blood: : 20 μ L

Qualitative determination of IgM and IgG antibodies to SARS-CoV-2 virus

Minimum 3 minutes to get results, generally within 10-15 minutes



Gold standard test for COVID-19
96 results within 2 hours
Simultaneous detection of the three genes of ORF1ab / E / N in 2019-nCoV

SARS-CoV-2 Nucleic Acid Detection Kit (Fluorescent PCR)

Step:1 Extract virus RNA from sample
Step:2 Perform RNA amplification

The reagent can be run on PCR analysers with 4 channels with FAM, ROX, Cy5 and HEX or VIC.
The Recommended analysers include: Thermo ABI 7500, ABI Q5,Roche Lightcycler 480II, Hongshi SLAN-96P,SLAN-96S).